

# ZOLL®



SeaSatCom

MARITIME CONNECTIVITY SOLUTIONS

## Are You Prepared for Sudden Cardiac Arrest at Sea?



Now more than ever, it's critical that you have an automated external defibrillator (AED) on board to provide immediate aid to a victim of sudden cardiac arrest (SCA).

It happens in a split second. A person collapses and the clock starts ticking. For every minute that passes without CPR, a victim's chance of survival decreases by 10%.

If bystanders wait for emergency services to arrive, cardiac arrest victims have just a 5% chance of survival. If those same bystanders immediately intervene with CPR and an AED, chance of survival increases to 60%.

SCA can happen to anyone at any time, anywhere — on land or sea. It doesn't discriminate and can affect otherwise healthy children and adults of all ages.



**350,000**

people suffer out-of-hospital cardiac arrests in Europe each year



**10%**

Decrease in survival for every minute of delayed defibrillation




**5%**

Chance of survival if a victim collapses and nothing is done until emergency services arrive

**FOR MORE INFORMATION: +33 4 89 73 01 32 or [contact@seasatcom.com](mailto:contact@seasatcom.com)**

SeaSatCom - 37 chemin du Puissanton - 06220 VALLAURIS- France - [www.seasatcom.com](http://www.seasatcom.com)

Sarl au capital de 39 390€ - RCS Antibes 482 017 787 - Siret 482 017 787 00037 - TVA intracom FR71 482 017 787



### ZOLL AED 3®


- Real CPR Help®
- Integrated paediatric rescue
- WiFi connectivity
- Long-life consumables



### AED Plus®

- Real CPR Help®
- Visual & audio prompts
- One-piece electrode
- Long-life consumables

Cardiac Arrest Survival Rate



## ZOLL AEDs offer real-time CPR feedback

Designed for unexpected heroes, all ZOLL® AEDs guide users through the process of performing high-quality CPR and will, if needed, deliver a potentially lifesaving shock to the heart.

Additionally, all ZOLL AEDs have an ingress protection (IP) rating of 55 and perform automatic self-testing, a series of regular tests that verify the AED, battery, and pads are fully functional.

## Accessories : Everything you need to keep your AED accessible and ready for use

### Zoll Medical Carrying Case Hard Shell For Zoll AED



### Zoll Medical Rescue Backpack For Zoll AED



### Lithium Manganese Dioxide battery pack



## SEASATCOM AED Program Management

- AED Program Management Onboard
  - The ZOLL AED 3 can be easily cloud-connected to ZOLL's program management system, PlusTrac, via Wi-Fi
  - Automatically reports status to PlusTrac
  - Failed self-test generates immediate e-mail
  - Alerts for needed replacement of battery or pads
  - Program managers can review all their AEDs any day, any time
  - PlusTrac
    - Complete AED program solution from ZOLL with Program Management Onboard
    - Readiness – Automated AED status and Trained Responder Tracking
    - Compliance – Consumable replacement and maintenance records
    - Seasatcom Support – Local customer with Post Use resources

### AED PROGRAM MANAGEMENT ONBOARD™



### Want to know more?

Please contact your ZOLL Medical Distributor :  
SEASATCOM - [contact@seasatcom.com](mailto:contact@seasatcom.com)

FOR MORE INFORMATION: +33 4 89 73 01 32 or [contact@seasatcom.com](mailto:contact@seasatcom.com)

SeaSatCom - 37 chemin du Puissanton - 06220 VALLAURIS- France - [www.seasatcom.com](http://www.seasatcom.com)

Sarl au capital de 39 390€ - RCS Antibes 482 017 787 - Siret 482 017 787 00037 - TVA intracom FR71 482 017 787